

ARE YOUR
Directional Fittings
On Thin Ice?



Small Footprint
w: 11" d: 8.5" h: 5"

Smaller than a typical test chamber set-up, **the MedRx AVANT Polar HIT™ Chamber** can be discreetly situated in any office setting

MedRx®

www.medrx-usa.com or Call Toll Free: 888-392-1234

AVANT

Polar HIT™

Directional Hearing Instrument Test Chamber



***Verify proper
functionality of directional
microphones prior to the
initial patient fitting***



- **Unique Rotating Test Pad
For Generating Directional
Polar Plots**
- **ANSI and IEC Hearing
Instrument Testing**
- **Simple, One-Cable USB
Connectivity**
- **Superior Acoustic Isolation
For Fast, Accurate Testing**

MedRx INC.®
www.medrx-usa.com

*1200 Starkey Rd. #105,
Largo, FL 33771
Toll Free 888-392-1234
727-584-9600 • Fax: 727-584-9602*

#0086 CE
AH-S-ADH093-2

MedRx's Authorized Representative in Europe
(Regulatory affairs only) Emergo Europe • Molenstraat 15
2513 BH The Hague • The Netherlands
Tel: (31) (0) 70 345-8570 • Fax: (31) (0) 70 346-7299